



Experience 3 weeks in Forschungzentrum

From 12.11.12 to 01.12.12

Mikheil Kelenjeridze

kelenjeridzemikheil@gmail.com

GTU

Faculty of Informatics and Control Systems

Advisor: Qetevan Kotetishvili

How we involved:

Begin in Tbilisi

Then was summer school in Batumi



• Then flight




Introducing day in FZJ



See IKP and COSY

Dr. Andro Kacharava explained us every necessary information









Were in **INM** and Ezequiel showed us their working labs

INM: Institute of Neuroscience and Medicine



\bigcirc C \bigcirc 1 1	C	• 4 •	C 11	1	(TXTX I 1)	`
Structural and	Hinciionai	organization	or the	nrain i		1
Structural and	i i anotionai	organization	OI tile	Oralli ((11 1111 1)	,


	Molecular organization of the brain	(INM-2)
--	-------------------------------------	---------

Cognitive Neuroscience (INM-3	3))
-------------------------------	----	---


Medical Imaging Physics	(INM-4)
-------------------------	---------

Nuclear Chemistry	(INM-5))
Tradition Silving	(11111	j


Ethics in the Neurosciences (INM-8)



DTI working principle



Isotropic diffusion



Anisotropic diffusion


By Einstein's equation:

$$\langle \Delta x^2(t) \rangle = 2Dt$$

Traveled distance in *t* time (what we measure with MRI)

Diffusion coefficient (what we infer from the MRI measurement)


Presenting information by ellipsoid




$$D = \begin{pmatrix} D_{xx} & D_{xy} & D_{xz} \\ D_{xy} & D_{yy} & D_{yz} \\ D_{xz} & D_{yz} & D_{zz} \end{pmatrix}$$

$$= \begin{pmatrix} \lambda_1 & 0 & 0 \\ 0 & \lambda_2 & 0 \\ 0 & 0 & \lambda_3 \end{pmatrix}$$


Presenting information by colors







Every ellipsoid color are related to its three main axis


Pixel color determines the ellipsoid's form





The last face of information











Phantoms for Studies of Anisotropic Diffusion in the Brain








Phantoms are tightly winded ribbons on plastic which will be placed in water and then analyzing water molecules diffusivity by MRI

What we were doing in INM

We made new experiment on phantom:

Winded ribbon on plastic tighter than on other phantoms



New phantoms:

Ribbon layers – 45

Ribbon rounds – 1880

Thickness – 19mm


After winding we cut some plastics for attaching put it in water box and then placed it in vacuum for removing air bubbles









MRI measurements










Getting information from MRI










How to determine fiber directions







Analyzing information









At the end of presentation

Thankful:

Prof. Dr. Hans Stroeher, for inviting us in FZJ and gave us chance for getting good experience

Dr. Andro Kacharava, for guiding us and introducing everything in FZJ

Prof. Ketevan Kotetishvili for leading us and giving necessary informations

Prof. Dr. Jon Shah, Dr. Farida Grinberg, Ezequiel Farrher and Dr. Helmut Soltner, for giving us opportunity involved in their practical activities

Thank you for your attention



Any questions?